

Optical clearing enabled by the Kramers-Kronig relation

Researchers in the Brongersma and Hong labs have developed a novel method for optical clearing of turbid biological tissues.

Studying how cells behave within the 3-dimensional context of the biological structures they make up is critical for fully understanding tissue and organ function. While optical imaging techniques are fundamental tools of biological and medical research, they are intrinsically limited to studying extremely thin sections of biological tissues due to how these tissues scatter light. The primary cause of this light scattering is the heterogeneous mix of refractive indices within biological tissues. Several techniques have been developed to correct this refractive index mismatch by replacing low-index water with high-index organics or vice versa; however, these techniques often involve toxic reagents and cannot be used in live tissues.

This new optical clearing method from the Brongersma and Hong labs leverages a distinct physical phenomenon compared to existing methods, features low-toxicity, and can even be used for optical clearing of live *in vivo* tissues and organs.




Fig. 1

Applications

- Optical/Tissue clearing of turbid biological tissues
 - Optical/Tissue clearing of live tissues
 - Optical/Tissue clearing of *in vivo* tissues

Advantages

- Low-toxicity
- applicable to live tissues
- applicable to *in vivo* tissues and organs

Innovators

- Guosong Hong
- Zihao Ou
- Nicholas Rommelfanger
- Shan Jiang
- Mark Brongersma
- Yi Shiou Duh
- Kenneth Brinson

- Su Zhao

Licensing Contact

Luis Mejia

Senior Licensing Manager, Physical Sciences

[Email](#)