

**Docket #:** S03-094

# **MRI - Broadband Multi-Coil Reconstruction**

Imaging time using PILS is reduced by using multiple coils with localized sensitivities with each coil having a separate demodulation channel thereby permitting parallel signal processing and image reconstruction. Images from the multiple coils are then combined to form an image with a larger field of view (FOV).

## **Innovators**

- Jin Hyung Lee
- John Pauly
- Dwight Nishimura

## **Licensing Contact**

### **Irit Gal**

Senior Licensing Manager

[Email](#)